



Single pack - P1001AD001

100-pack - P1001AD011



## SalivaPOD<sup>®1</sup>

### Straw-based Saliva Collection for Optimal Assay Performance

WHEN SAMPLING TIME, CAPACITY AND PATIENT COMPLIANCE COUNTS

#### Intended Purpose (Area of Use)

“The SalivaPOD<sup>®</sup> is a single-use device intended for use by a lay person for the collection and short-term storage of oral fluid specimens. The lay person may perform, or the healthcare professional performs, the transfer of the specimen into a transport culture medium, a buffer solution or one or more diagnostic reagents for further transport to a testing laboratory, or for immediate analysis and diagnostic use.”



ConceptoMed

# SalivaPOD<sup>®1</sup> - Collector and Injector

## Description, information and benefits

- Premium STRAW-based saliva sampling device delivering low viscosity samples
- Easy-to-use for the home setting
- Standardizes saliva collection for delivery of optimal assay performance in lab processing
  - ▶ Low variance in saliva-to-preservative ratio due to controlled and precise titration of optimal volume of saliva
  - ▶ The closed system transfer to preservative tube eliminates risk of cross contamination and false test results
- Minimizes work load for healthcare professionals
- Substantial improvement of patient convenience, compared to funnel based systems

## Structural design elements

- The SalivaPOD™ Mouthpiece supports pipetting with robot by allowing low viscosity saliva to enter the volume chamber
- 8 arrows and '360° Fill-to-Line mark' collects ca. 1 ml saliva, with air bubbles above line
- Over-flow threshold guarantees maximum 1.4 ml saliva may be injected into 16 mm vial
- Every label prepared with printed TWIN Unique QR code for integration with eHealth systems
- Easy-to-open IFU Package Insert with 13 languages
- Available interactive eIFU SalivaCODE™ App with optional oral reading of text
- Optional peel-off label on the IFU paper insert (single version only) for manual written marking of sample

## Material specifications - technical info

Volume chamber:	☐ Transparent Polypropylene <sup>2</sup>
Mouthpiece:	☐ ABS <sup>2</sup>
Protection cap:	■ Polycarbonate/ABS <sup>2</sup>
Bottom screw cap:	■ Polypropylene <sup>2</sup>


Mouthpiece top part :	Female luer slip connector
Bottom Screw Threads:	Fits 16mm tubes from Copan, BD ++
Fill-to-Line (bubbles above line):	Approx. 1 ml saliva injected into vial
Over-flow threshold:	Max. 1.4 ml saliva injected into vial

## Primary device packaging

Material	High quality foil and Flexpeel® Tyvek 2FS CR27 - Strong material supporting proper storage and use
Blister	Blister stays open and does not return to the initial configuration after opening (no memory).

## IFU paper insert

About	Paper Package insert with alternative sticker label 32 mm x approx 70 mm x 7 mm
Languages	Norwegian, Swedish, Danish, Finnish, English, Spanish, German, French, Italian, Czech, Polish, Somali, Arabic
Content	Step 1. Saliva sampling procedure Step 2. Transfer of saliva into a compatible laboratory vial
eIFU	PRINT READY IFU (A4 size) can also be downloaded from website <a href="http://ifu.salivapod.com">ifu.salivapod.com</a> for distribution SalivaCODE App can be used as eIFU for the saliva collection procedure eIFU and downloadable IFUs can be expanded with selected languages



# SalivaPOD<sup>®1</sup> - Content and Packaging

## Single Pack

### Secondary packaging

- Material** Cardboard material. Outer dimensions:  
39 mm x 38 mm x 145 mm
- Label** IFU and secondary label on the side of the secondary box
- Content in box** Primary packing with SalivaPOD™  
IFU paper insert, 13 languages

### Outer packaging (tertiary)

- Material** Cardboard material. Outer dimensions:  
300 mm x 207 mm x 217 mm
- Label** Tertiary case label on the side of the tertiary box
- Content in box** 50 pcs SalivaPOD™ secondary pack, single-packs per case

## 100-pack

### Secondary packaging

- Material** Closed plastic bag with 100 pcs SalivaPOD in soft blister inside cardboard tertiary packaging (cardboard box)
- IFU** 1 pc IFU insert outside the secondary pack, for illustration  
Bundle of 100 pcs "Quick-Guide for collection procedure" incl. by default (Norwegian + English)

### Outer packaging (tertiary)

- Material** Corrugated 3mm cardboard material. Outer dimensions:  
300 mm x 207 mm x 217 mm
- Content in box** 1 sec.pack. bag of 100 pcs. SalivaPOD devices in primary packaging  
1 single IFU for reference (outside sec.pack)  
Bundle of 100 pcs "Quick-Guide for collection procedure" incl. by default (Norwegian + English)




## User steps



## SalivaPOD™ screw-fits onto vials in broad use today

SalivaPOD™ fits seamlessly into existing diagnostics hardware compatible vials (Ø 16mm tubes from e.g. COPAN, BD ++)



# SalivaPOD<sup>®1</sup> - SKUs, Ordering details



	Product description	Art.no	Quantity (cases x pcs.)	SalivaPODs per half-pallet	W cm	D cm	H cm	Weight (kg.)
Half-pallet	SalivaPOD™ 1 blister pack 1 IFU	P1001AD001	60 cases x 50	3 000	80	120	105	105
Half-pallet	SalivaPOD™ 100 blister pack 1 IFU 100 simple IFUs	P1001AD011	60 cases x 100	6 000	80	120	105	155

## Pallets (EUR/EPAL)

Full pallet: 2 half-pallets, stacked  
Size: 80cm x 120cm x (105cm x 2)

## Manufacturing - country of origin

Designed in Norway by ConceptoMed AS  
Manufactured and assembled in Sweden

# SalivaCODE™ - App service - Optional LIMS API

SALIVA  
CODE



## eIFU SalivaCODE™ is available

- Works with SalivaPOD™ QR codes.
- Interactive App/eIFU guides the user, ensuring optimal collection.
- Printed text/images/videos and optional oral reading of text
- Includes 13 languages today
- API available for integration of QR code usage to link saliva specimen with Laboratory Information Management Systems



SalivaPOD<sup>®</sup> is widely protected by granted and pending patents in USA, Canada, Europe, Asia

ConceptoMed

1. CE marked, Class 1, single use, non-sterile medical device, according to MDR. FDA listed as exempt device.  
2. Material Data Sheets on file - all materials are Medical Grade.

FDA listed, reg.no: REN# 3013576743, ON# 10054844

